

For trace contaminants in drinking water and wastewater

On-Line SPE LC/MS/MS Applications for Environmental Analysis

Xavier Maes
Agilent Technologies

The Measure of Confidence

Agenda Online SPE Solution

- Solid phase extraction
- Setup of automated online SPE
- Example 1: Polycyclic Aromatic Hydrocarbons
- Example 2: Pesticide
- Conclusion

Agilent Technologies

Sample Preparation

Goals for liquid Chromatography

Matrix removal (dissolved and non-dissolved constituents)

- Protects analytical columns
- Increases method selectivity and robustness
- Improves quantitation

Analyte enrichment

- Improves trace analysis performance
- Makes analysis feasible, if original concentration is below limit of detection (LOD)

Sample solvent adaptation (if limited compatibility with LC method)

- Makes analysis feasible
- Improves resolution, LOD, robustness

Agilent Technologies

Sample Preparation

Techniques

Most common techniques

- Solid phase extraction (SPE)
- Liquid-liquid extraction (LLE)
- Precipitation (e.g., proteins) and centrifugation
- Membrane filtration, ultra centrifugation

Agilent Technologies

Disadvantages of offline SPE

Often performed manually, there are several disadvantages with this approach.

- Time-consuming, labor and cost intensive, and limits productivity
- Exposure to hazardous or infectious matrices (such as biofluids) raises safety concerns
- Batch to batch reproducibility of the recovery of the analyte can vary

Agilent Technologies

Typical sample preparation for water

- For water screening and quantitation in the 10 ppt range usually requires solid phase extraction (SPE) to pre-concentrate the target compounds in the sample solution

10 ng/L
1000 mL water sample

1 ug/L
1 mL organic extract

10 pg o.c.
Injection volume: 10 µL


Agilent Technologies

Online SPE


Online Solid Phase Extraction

Automating this SPE process removes the limitations of the manual SPE process:

- Enables direct injection of untreated samples
- Automates sample cleanup and analyte enrichment
- Automation allows for a faster overall workflow
- higher batch to batch reproducibility
- Samples can be run unattended.



Online SPE Process




Valve Position 1


- Load sample 1 to SPE 1

Features

- Online SPE with 2 x SPE cartridges
- Injection volume up to 1800 µl (2 x 900 µl)




Online SPE Process




Valve Position 2

- Elute sample 1 from SPE 1 to LC column
- Load sample 2 to SPE 2




Application 1

Sensitive Determination of Polycyclic Aromatic Hydrocarbons in Tap Water by Online Solid Phase Extraction and UHPLC



Polycyclic Aromatic Hydrocarbons

- Formed by the incomplete burning and combustion of hydrocarbons.
- Adsorb onto flying particles which can form as rain drops and snow.
- Leaching from linings of water storage tanks and distribution lines.
- 7 PAHs are classified as probable human carcinogens by the Environmental Protection Agency (US EPA)




Polycyclic Aromatic Hydrocarbons Analysis

EPA / EU Regulations

According to EPA (method 550, 550.1 and 610) and the European Union legislation, 16 parent PAHs should be monitored.


1 Naphthalene	5 Phenanthrene	9 Benzo(a)anthracene	13 Benzo(a)pyrene
2 Acenaphthylene	6 Anthracene	10 Chrysene	14 Dibenzo(a,h)anthracene
3 Acenaphthene	7 Fluoranthene	11 Benzo(b)fluoranthene	15 Benzo(g,h,i)perylene
4 Fluorene	8 Pyrene	12 Benzo(k)fluoranthene	16 Indeno(1,2,3-cd)pyrene

- based on offline SPE or LLE
- require **high sample volumes (1 L)**
- are **labor intensive and time-consuming**
- have **run times greater than 60 minutes**



Chromatographic Conditions



Parameter	Value
Analytical column	Agilent ZORBAX Eclipse PAH 4.8 × 150 mm, 3.5 µm (p/n 959991-918)
Column temperature	25 °C
Gradient	A) Water, B) Acetonitrile 0 minutes, 40 %B at 1.3 mL/min 5 minutes, 40 %B 20 minutes, 100 %B 30 minutes, 100 %B 30.1 minutes, 40 %B Posttime 6 minutes
Injection volume	1,800 µL (online SPE), needle wash 3 seconds with acetonitrile
Autosampler temperature	18 °C
Online SPE enrichment column	Polaris 5 C18-A, 3.0 × 50 mm (p/n A290090D930)
Online SPE enrichment method	A) Water, B) Acetonitrile 0 minutes, 20 %B at 1.3 mL/min 3 minutes, 20 %B 4 minutes at 0 mL/min 5 minutes, 100 %B at 1 mL/min 9 minutes, 100 %B 9.1 minutes, 20 %B at 1 mL/min 15 minutes, 20 %B at 0 mL/min 15.1 minutes, 20 %B at 0 mL/min
Diode array detector	230 nm, bandwidth 4 nm, reference 400 nm, reference bandwidth 100 nm, 10 Hz
Fluorescence detector	Multisignal acquisition, set at λ _{ex} = 280 nm and λ _{em} = 350 nm (FLD A), 330 nm (FLD B), 440 nm (FLD C), 500 nm (FLD D), 11.51 Hz, PMT 13



Chromatographic Results

RT, RT & Area RSDs, Detection

Elution order	Compound	Retention time (min)	Detection (nm)	RT RSD (%) (n = 7)	Area RSD (%) (n = 7)
1	Naphthalene	10.6	UV 230/Ref 400	0.028	3.8
2	Acenaphthylene	11.4	UV 230/Ref 400	0.025	0.9
3	Acenaphthene	12.5	Ex 260/Em 330	0.029	1.0
4	Fluorene	12.8	Ex 260/Em 330	0.031	0.5
5	Phenanthrene	13.5	Ex 260/Em 350	0.035	0.9
6	Anthracene	14.3	Ex 260/Em 440	0.040	0.6
7	Fluoranthene	15.1	Ex 260/Em 440	0.050	0.3
8	Pyrene	15.7	Ex 260/Em 440	0.050	6.0
9	Benzo(a)anthracene	17.7	Ex 260/Em 440	0.070	2.3
10	Chrysene	18.2	Ex 260/Em 350	0.071	0.8
11	Benzo(b)fluoranthene	19.7	Ex 260/Em 440	0.080	1.3
12	Benzo(k)fluoranthene	20.5	Ex 260/Em 440	0.090	1.1
13	Benzo(a)pyrene	22.1	Ex 260/Em 440	0.080	2.5
14	Dibenzo(a,h)anthracene	22.2	Ex 260/Em 440	0.011	2.0
15	Benzo(g,h,i)perylene	22.8	Ex 260/Em 440	0.013	1.4
16	Indeno(1,2,3-cd)pyrene	23.5	Ex 260/Em 500	0.015	3.9





Application 2

Quantitative LC/MS/MS Analysis of 22 pesticides including Dichlorophenols in Water Using on-line SPE Enrichment.

Results requirement

Name	LOQ µg/L	CV %	Recovery %
2,4-D	0.010	5	98-102
MCPA	0.010	5	98-102
Mechlorprop	0.010	5	98-102
Dinoseb	0.010	5	98-102
2,6-dichlorophenol *	0.010	5	98-102
Terbutylazin	0.010	5	98-102
Isoproturon	0.010	5	98-102
Altrazin	0.010	5	98-102
Bentazone	0.010	5	98-102
Pendimethalin	0.010	5	98-102
Cymazin	0.010	5	98-102
2-Hydroxy-atrazin	0.010	5	98-102
Desethylatrazin	0.010	5	98-102
Hexazinon	0.010	5	98-102
DNOC	0.010	5	98-102
Dichlorprop	0.010	5	98-102
Desisopropylatrazin	0.010	5	98-102
2,6-dichlorbenzamid *	0.010	5	98-102
Dimethoat	0.010	5	98-102
Metamitron	0.010	5	98-102
2,4-dichlorophenol *	0.010	5	98-102
Simazin	0.010	5	98-102

Agilent Technologies



- ### Summary and Conclusions
- Online SPE/LC-QQQ simple, fast and reliable analysis at concentration of ng/L range.
 - fully controlled and automate
 - Good recovery reproducibilities
 - time and labor savings compared to off line SPE
- Agilent Technologies

Thank you for your Attention

Agilent Technologies